















Choose with confidence

Motiva Anatomical TrueFixation™

 (cm)	 (cm)	 (cm)	 (cm)	V (cc)	 (cm)	 (cm)	 (cm)	 (cm)	V (cc)
FF FULL HEIGHT FULL PROJECTION 					FC FULL HEIGHT CORSÉ PROJECTION 				
11.0	11.5	4.4	6.8	325	11.0	11.5	5.4	7.7	400
11.5	12.0	4.6	7.1	370	N/A	N/A	N/A	N/A	N/A
12.0	12.5	4.8	7.4	420	12.0	12.5	5.8	8.3	505
12.5	13.0	5.0	7.7	475	12.5	13.0	6.0	8.6	560
13.0	13.5	5.2	8.0	535	N/A	N/A	N/A	N/A	N/A
13.5	14.0	5.4	8.3	600	N/A	N/A	N/A	N/A	N/A
MF MEDIUM HEIGHT FULL PROJECTION 					MC MEDIUM HEIGHT CORSÉ PROJECTION 				
11.0	10.0	4.4	6.4	280	11.0	10.0	5.4	7.3	340
11.5	10.5	4.6	6.7	320	11.5	10.5	5.6	7.6	390
12.0	11.0	4.8	7.0	365	12.0	11.0	5.8	7.9	440
12.5	11.5	5.0	7.3	415	12.5	11.5	6.0	8.2	495
13.0	12.0	5.2	7.6	465	13.0	12.0	6.2	8.5	555
13.5	12.5	5.4	7.9	525	N/A	N/A	N/A	N/A	N/A
LF LOW HEIGHT FULL PROJECTION 					LC LOW HEIGHT CORSÉ PROJECTION 				
11.0	9.0	4.4	6.2	250	11.0	9.0	5.4	7.1	305
11.5	9.5	4.6	6.5	285	11.5	9.5	5.6	7.4	350
12.0	10.0	4.8	6.8	330	12.0	10.0	5.8	7.7	400
12.5	10.5	5.0	7.1	375	12.5	10.5	6.0	8.0	450
13.0	11.0	5.2	7.4	425	N/A	N/A	N/A	N/A	N/A
13.5	11.5	5.4	7.7	480	N/A	N/A	N/A	N/A	N/A

A=Base D=Height B=Projection C=Inferior Arc Distance V=Volume

Anatomical TrueFixation™

Full	Corsé	SilkSurface®	ProgressiveGel PLUS™	BluSeal®	Single Pack	TwinPack™	Sizers	Q Inside	5Y Extended Warranty*
✓	✓	✓	✓	✓	✓	✓	✗	✓	✓

*Optional

